



**INSTRUMENT
TABLE**

Specification

Minimum Height	560mm
Maximum Height	860mm
Allowable Load	80 kg
Power	AC100-240V 50/60 Hz
Size of Base	604mm x 550 mm
Size of Table top	600mm x 400 mm
Battery Powered	N/A
Weight	20 kg

Description

The LightBox Instrument Table features motorized height adjustment, a new ergonomic design, and smooth caster movements with locks. It is available with a single or double tabletop, and can support loads up to 110 lbs.

Key Features


New ergonomic design

Motorized height adjustment

Smooth caster movements with locking mechanisms

Accepts loads up to 110 lbs. (50kg)

Single or double tabletop available



Add: Warehouse L1-05, Sheikh Rashid Bin Saeed Al Maktoum Street,
City: Ajman, UNITED ARAB EMIRATES
Tel +971-588051849